

APRIL peptide

Cat. No. GTX00459

Application	Neutralizing/Inhibition
Species	Human

Package
50 µg

APPLICATION

Application Note

The peptide is used for blocking the antibody activity of anti-APRIL (cat. No GTX31268). It usually blocks the antibody activity completely in Western blot by incubating the peptide with equal volume of antibody for 30 min at 37°C.

PROPERTIES

Form	Liquid
Buffer	PBS, 0.1% BSA
Preservative	0.02% Sodium azide
Storage	Store the peptide at -20°C, stable for one year.
Concentration	0.20 mg/ml (Please refer to the vial label for the specific concentration.)
Conjugation	Unconjugated
Note	For <i>In vitro</i> laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).