

F. t. holarctica FTL0187 antibody [PAb (196)]

Cat. No. GTX00495

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA
Reactivity	Francisella tularensis ssp. holarctica

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Product Note

The antibody does not react with Francisella tularensis subsp. tularensis homolog.

Properties

Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant FLT0187 (cyclohexadienyl dehydratase) protein produced in E. coli
Purification	Purified by affinity chromatography From serum
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX00495 WB Image

WB analysis of various samples using GTX00495 *F. t. holarctica* FTL0187 antibody [PAb (196)].

St : Mw marker

FSC200-CL : 10 µg of FSC200 cell lysate

SCHUS4-CL : 10 µg of SCHU S4 cell lysate

FSC200-MF : 10 µg of FSC200 membrane fraction

SCHUS4-CL : 10 µg of SCHU S4 membrane fraction

K1 : 50 ng of FTL0187 recombinant protein

K2 : 125 ng of FTL0187 recombinant protein



For full product information, images and publications, please visit our [website](https://www.genetex.com).