

# F. t. tularensis FTT0571 antibody [PAb (181)]

**Cat. No. GTX00501**

|                     |                                        |
|---------------------|----------------------------------------|
| <b>Host</b>         | Rabbit                                 |
| <b>Clonality</b>    | Polyclonal                             |
| <b>Isotype</b>      | IgG                                    |
| <b>Applications</b> | WB, ELISA                              |
| <b>Reactivity</b>   | Francisella tularensis ssp. tularensis |

**Package**  
100 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | Assay dependent      |
| ELISA              | Assay dependent      |

Not tested in other applications.

### Product Note

The antibody does not react with Francisella tularensis subsp. holarctica homolog.

## Properties

|                      |                                                                                                        |
|----------------------|--------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                 |
| <b>Buffer</b>        | PBS                                                                                                    |
| <b>Preservative</b>  | 15mM Sodium azide                                                                                      |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. |
| <b>Concentration</b> | 1 mg/ml (Please refer to the vial label for the specific concentration.)                               |
| <b>Immunogen</b>     | Recombinant FTT0571 protein produced in E. coli                                                        |
| <b>Purification</b>  | Purified by affinity chromatography<br>From serum                                                      |
| <b>Conjugation</b>   | Unconjugated                                                                                           |

### Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX00501 WB Image

WB analysis of various samples using GTX00501 *F. t. tularensis* FTT0571 antibody [PAb (181)].

St : Mw marker

FSC200-CL : 10 µg of FSC200 cell lysate

SCHUS4-CL : 10 µg of SCHU S4 cell lysate

FSC200-MF : 10 µg of FSC200 membrane fraction

SCHUS4-CL : 10 µg of SCHU S4 membrane fraction

K1 : 50 ng of FTT0571 recombinant protein

K2 : 125 ng of FTT0571 recombinant protein



For full product information, images and publications, please visit our [website](https://www.genetex.com).