

Mouse Anti-Human kappa light chain antibody [TB28-2] (APC-Cy7)

Cat. No. GTX00506-15

Host	Mouse
Clonality	Monoclonal
Isotype	lgG1
Applications	FCM
Reactivity	Human

Package 100 test

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	4 μ l reagent / 100 μ l of whole blood or 10 6 cells in a suspension
Not tested in other applications.	
Product Note	This antibody recognizes human Ig kappa light chains of both secreted and cell surface immunoglobulin. It detects also free kappa light chains.
Properties	
Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Human IgG-kappa myeloma protein
Purification	Purified by size-exclusion chromatography
Conjugation	Allophycocyanin-Cyanine7 (APC-Cy7) Wavelength
Note	For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Cy and CyDye are trademarks of GE Healthcare.



For full product information, images and publications, please visit our website.

Date 2025 / 10 / 21 Page 1 of 1