

# KIR2DL5A antibody [UP-R1]

**Cat. No. GTX00515**

|                     |             |
|---------------------|-------------|
| <b>Host</b>         | Mouse       |
| <b>Clonality</b>    | Monoclonal  |
| <b>Isotype</b>      | IgG1        |
| <b>Applications</b> | ICC/IF, FCM |
| <b>Reactivity</b>   | Human       |

**Package**  
100 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| ICC/IF             | Assay dependent      |
| FCM                | Assay dependent      |

Not tested in other applications.

## Properties

|                      |                                                                                                        |
|----------------------|--------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                 |
| <b>Buffer</b>        | PBS                                                                                                    |
| <b>Preservative</b>  | 15mM Sodium azide                                                                                      |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. |
| <b>Concentration</b> | 1 mg/ml (Please refer to the vial label for the specific concentration.)                               |
| <b>Immunogen</b>     | Human CD158f-Ig fusion protein                                                                         |
| <b>Purification</b>  | Protein A purified                                                                                     |
| <b>Conjugation</b>   | Unconjugated                                                                                           |


### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX00515 FCM Image**

FACS analysis of human peripheral blood using GTX00515 KIR2DL5A antibody [UP-R1].



For full product information, images and publications, please visit our [website](#).