

# GRB2 antibody [PAb (554)]

**Cat. No. GTX00529**

|                     |            |
|---------------------|------------|
| <b>Host</b>         | Rabbit     |
| <b>Clonality</b>    | Polyclonal |
| <b>Isotype</b>      | IgG        |
| <b>Applications</b> | WB         |
| <b>Reactivity</b>   | Human      |

**Package**  
100 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | Assay dependent      |

Not tested in other applications.

**Calculated MW** 25 kDa. ( [Note](#) )

## Properties

|                      |                                                                                                        |
|----------------------|--------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                 |
| <b>Buffer</b>        | PBS                                                                                                    |
| <b>Preservative</b>  | 15mM Sodium azide                                                                                      |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. |
| <b>Concentration</b> | 1 mg/ml (Please refer to the vial label for the specific concentration.)                               |
| <b>Immunogen</b>     | Synthetic peptide corresponding to C-terminus of human Grb2 (amino acids 198-217).                     |
| <b>Purification</b>  | Purified by affinity chromatography<br>From serum                                                      |
| <b>Conjugation</b>   | Unconjugated                                                                                           |

### Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX00529 WB Image

WB analysis of Ramos cell lysate using GTX00529 GRB2 antibody [PAb (554)].



For full product information, images and publications, please visit our [website](#).