

ZAP70 antibody [ZAP-03]

Cat. No. GTX00530

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, FCM
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
FCM	Assay dependent

Not tested in other applications.

Calculated MW 70 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Bacterially expressed fusion protein representing C-terminal part (160 amino acids) of human ZAP70 with histidine tag
Purification	Protein A purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTx00530 WB Image

WB analysis of HPB-ALL cell lysate using GTx00530 ZAP70 antibody [ZAP-03].

Lane 1, 2, 3, 4 : GTx00530 (Dilution : 4, 2, 1, 0.5 µg/ml)

Lane 5 : anti-ZAP-70 comparative antibody

Lane 6 : mouse IgG1 isotype control


GTx00530 WB Image

WB analysis of HPB-ALL (Lane 1, 3) and Ramos cell (Lane 2, 4) lysates using GTx00530 ZAP70 antibody [ZAP-03].

Lane 1, 2 : GTx00530

Lane 3, 4 : mouse IgG1 isotype control

Dilution : 0.5µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).