

# SSEA-1 / CD15 antibody [MMA] (PE)

## Cat. No. GTX00545-08

Host	Mouse
Clonality	Monoclonal
Isotype	IgM
Applications	FCM
Reactivity	Human

Package 100 test

## Applications

#### **Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution Recom	mended dilution
FCM 10 µl re	eagent / 100 μl of whole blood or 10 <sup>6</sup> cells in a suspension

Not tested in other applications.


Properties	
Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	U937 histiocytic lymphoma cells
Purification	Purified by size-exclusion chromatography
Conjugation	Phycoerythrin (PE) Wavelength
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 12 / 12 Page 1 of 2

### DATA IMAGES



#### GTX00545-08 FCM Image

FACS analysis of human peripheral blood using GTX00545-08 CD15 antibody [MMA] (PE).



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 12 / 12 Page 2 of 2