

CD4 antibody [OX-35] (APC)

Cat. No. GTX00567-07

Host	Mouse
Clonality	Monoclonal
Isotype	lgG2a
Applications	FCM
Reactivity	Rat

Package 100 μg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

	•
Suggested dilution	Recommended dilution
FCM	Assay dependent
Not tested in other applications.	
Product Note	This antibody reacts with an extracellular epitope of rat CD4 transmembrane glycoprotein (55 kDa).
Properties	
Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	MLR generated rat Th cells
Purification	Purified by size-exclusion chromatography
Conjugation	Allophycocyanin (APC) Wavelength
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our <u>website</u>.

determine the structure or sequence of the product.

Date 2025 / 11 / 23 Page 1 of 1

€ 886-3-6208988 📻 886-3-6208989 🐷 infoasia@genetex.com

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to