

STYK1 antibody

Cat. No. GTX00635

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent

Not tested in other applications.

Calculated MW 48 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 2% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5-1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human STYK1
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES



GTX00635 WB Image

WB analysis of HepG2 cell lysate using GTX00635 STYK1 antibody.


Dilution : 0.2-1 µg/ml



GTX00635 IHC-P Image

IHC-P analysis of human colon tissue using GTX00635 STYK1 antibody.

Dilution : 1:100



GTX00635 IHC-P Image

IHC-P analysis of human brain, cerebellum tissue using GTX00635 STYK1 antibody.

Dilution : 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).