

EID1 antibody [2]

Cat. No. GTX00708

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ICC/IF
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	~1 µg/ml
ICC/IF	~1 µg/ml

Not tested in other applications.

Calculated MW 21 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Filter-sterilized PBS, 50% Glycerol
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide containing amino acids 1-19 of human EID1 protein
Purification	Purified IgG From tissue culture supernatant
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX00708 ICC/IF Image

ICC/IF analysis of HeLa cells using GTX00708 EID1 antibody [clone 2].


Red : Primary antibody

Green : DAPI

Dilution : 1:1000

Fixation : 4% PFA, overnight

Permeabilization : 0.25% Triton X-100 in PBS for 10 min



GTX00708 WB Image

WB analysis of MCF-7 mock (Lane 1) or overexpressing EID1 (Lane 2) using GTX00708 EID1 antibody [clone 2].

The EID1 protein was identified as a 21kDa protein band as shown by an arrow.



For full product information, images and publications, please visit our [website](https://www.genetex.com).