

IGSF8 antibody

Cat. No. GTX00745

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000
ICC/IF	1:100
IHC-P	1:100

Not tested in other applications.

Calculated MW 65 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full-length mouse IGSF8 with Flag tag
Purification	Purified IgG
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX00745 WB Image


WB analysis of various mouse tissue samples of Igff8-targeted mice using GTX00745 IGSF8 antibody.

(F) : floxed alleles

(-) : knock-out alleles

Dilution : 1:500

Loading : 30 µg




GTX00745 ICC/IF Image

ICC/IF analysis of NIH-3T3 cells using GTX00745 IGSF8 antibody.

Dilution : 1:100

Fixation : 4% PFA (v/v)




GTX00745 WB Image

WB analysis of NIH-3T3 whole cell lysate using GTX00745 IGSF8 antibody.

Dilution : 1:1000

Loading : 40 µg

SDS-PAGE : 12.5%



GTX00745 IHC-P Image

IHC-P analysis of mouse lung tissue section using GTX00745 IGSF8 antibody.

Dilution : 1:100

Antigen retrieval : Proteolytic-Induced Epitope Retrieval



For full product information, images and publications, please visit our [website](https://www.genetex.com).