

PDILT antibody

Cat. No. GTX00747

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, IHC-Fr, IP
Reactivity	Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
IHC-P	1:100-1:500
IHC-Fr	1:100-1:500
IP	1:100

Not tested in other applications.

Calculated MW 68 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Serum
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Two synthetic peptides corresponding to the C-terminal regions of mouse PDILT, C+IRKPEEPERRKETA (550-563) and C+QPKEQPKPERKLEV (571-584), respectively, conjugated with KLH
Purification	Unpurified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX00747 IHC-Fr Image

IHC-Fr analysis of mouse testicular tissue section (5 μ m) using GTX00747 PDILT antibody.

Dilution : 1:100

Fixation : 4% PFA in PBS


GTX00747 WB Image


WB analysis of various mouse tissue samples using GTX00747 PDILT antibody.

This protein is testis-specific.

Dilution : 1:1000

Loading : 30 μ g for Tissue proteins & 10 μ g for sperm protein

This protein is testis-specific.


GTX00747 IP Image

IP analysis of wild-type (+/+) and Pdilt (-/-) mouse testis samples using GTX00747 PDILT antibody.

Lane 1 : Input


Lane 2 : Wild-type

Lane 3 : Pdilt (-/-)

Total protein amount : 100 μ g

IP antibody : 1:100

IB antibody : 1:1000


GTX00747 IHC-P Image

IHC-P analysis of mouse testis tissue section using GTX00747 PDILT antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).