

Human Parathyroid (Normal) tissue slides

Cat. No. GTX00762

Application	IHC-P
Species	Human

Package
5 slide

PRODUCT

Summary

Tissue sections are 5 µm in thickness, and mounted on positively charged glass slides. Tissue was fixed in formalin immediately after excision, and embedded in paraffin. All GeneTex's paraffin embedded tissue products including tissue array, tissue sections, and tissue panels can be used for both immunohistochemistry and in-situ hybridization. These products are ideal for rapid cellular localization of RNA and protein expression. At least one of the tissue slides from each lot was stained with H & E to ensure the quality.

Tissue type Normal

APPLICATION

Application Note

Please bake slides at 60°C for 30 minutes before use.

PROPERTIES

Storage Cool dry at room temperature. 4°C for long term storage.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).