

AU5 tag antibody (Agarose)

Cat. No. GTX00767

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Application	IP
Reactivity	Species independent

Package
1 mg

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IP	15-25 µl antibody per 0.1-1 mg of protein extract or lysate

Note : Dilute in PBS or medium that is identical to that used in the assay system.

Not tested in other applications.

PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	AU5 (TDFYLK) conjugated to KLH
Purification	Purified by antigen-affinity chromatography
Conjugation	Agarose Antibody is coupled to agarose beads with cyanogen bromide.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).