

Mineralocorticoid Receptor antibody

Cat. No. GTX00779

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.1-0.5µg/ml

Not tested in other applications.

Calculated MW 107 kDa. ([Note](#))

Product Note Superfamily members of NR3C2 are not reactive to this antibody.

Properties

Form	Liquid
Buffer	5mg BSA, 0.9mg NaCl, 0.2mg Na ₂ HPO ₄
Preservative	0.05mg Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	This antibody was raised against a synthetic peptide corresponding to a sequence at the C-terminus of human NR3C2 (950-984aa HALKVEFPAMLVEIISDQLPKVESGNAKPLYFHRK), different from the related mouse sequence by one amino acid, and from the related rat sequence by two amino acids.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx00779 WB Image


WB analysis of various samples using GTx00779 Mineralocorticoid Receptor antibody.

Lane 1 : Mouse Brain Tissue Lysate

Lane 2 : Rat Brain Tissue Lysate

Dilution : 0.5 µg/mL

Loading : 50µg



GTx00779 WB Image

WB analysis of various samples using GTx00779 Mineralocorticoid Receptor antibody.

Lane 1 : Rat Kidney Tissue Lysate

Lane 2 : Mouse Kidney Tissue Lysate

Lane 3 : HeLa Whole Cell Lysate

Lane 4 : A431 Whole Cell Lysate

Dilution : 0.5 µg/mL

Loading : 50µg



For full product information, images and publications, please visit our [website](https://www.genetex.com).