

Histone H3 antibody [1B1-B2]

Cat. No. GTX01093

Host	Mouse
Clonality	Monoclonal
Isotype	IgG3
Applications	WB, IP, ChIP assay
Reactivity	Human, Mouse, Yeast

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
IP	Assay dependent
ChIP assay	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	C-terminal peptide of human histone H3
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX01093 ChIP assay Image

ChIP analysis of sample using GTX01093 Histone H3 antibody [1B1-B2].


GTX01093 WB Image

WB analysis of HeLa histone lysates using GTX01093 Histone H3 antibody [1B1-B2].


GTX01093 ChIP assay Image

ChIP analysis of sample using GTX01093 Histone H3 antibody [1B1-B2].



For full product information, images and publications, please visit our [website](https://www.genetex.com).