

Mouse CDC2 protein, GST and His tag

Cat. No. GTX01208-pro

Species Mouse

Package
50 µg

APPLICATION

Calculated MW 34 kDa. ([Note](#))

Observed MW 57 kDa.

PROPERTIES

Form Liquid

Buffer 20mM Tris (pH8.0), 500mM NaCl, 10% Glycerol

Preservative No preservative

Storage Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Concentration Batch dependent (Please refer to the vial label for the specific concentration.)

Region/Sequence N-terminal GST and poly His tagged; Met1-Met297 of Mouse CDC2 protein (P11440)

Expression System Baculovirus (Insect cells)

Purity > 85% by SDS-PAGE.

Endotoxin < 1.0 EU/µg (determined by LAL method)

Conjugation Unconjugated


Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX01208-pro Image

SDS-PAGE of 5 μ g GTX01208-pro Mouse CDC2 protein, GST and His tag.



For full product information, images and publications, please visit our [website](#).