

Mouse MMP9 protein

Cat. No. GTX01246-pro

Species	Mouse
---------	-------

References (1)

Package

20 µg

Applications

Calculated MW	81 kDa. (Note)
---------------	----------------------------------

Observed MW (kDa)	100 kDa.
-------------------	----------

Properties

Form	Lyophilized powder
------	--------------------

Buffer	Batch dependent (please contact us for details)
--------	---

Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Avoid freeze-thaw cycles.
---------	---

Region/Sequence	No tagged; Met1-Pro730 of Mouse MMP9 protein (P41245-1)
-----------------	---

Expression System	HEK293 cells
-------------------	--------------

Purity	> 95% by SDS-PAGE.
--------	--------------------

Endotoxin	< 1.0 EU/µg (determined by LAL method)
-----------	--


Conjugation	Unconjugated
-------------	--------------

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX01246-pro Image

SDS-PAGE of 5 µg GTX01246-pro Mouse MMP9 protein.



For full product information, images and publications, please visit our [website](#).