

ICOS Ligand antibody [HK5.3]

Cat. No. GTX01490

Host	Rat
Clonality	Monoclonal
Isotype	IgG2a
Applications	FCM, Neutralizing /Inhibition
Reactivity	Mouse

References (3)

Package

100 µg

PRODUCT

Summary

The HK5.3 antibody reacts with mouse CD275, an approximately 40 kDa protein that is a member of the Ig superfamily. CD275 is also known as B7-H2, ICOS ligand, and B7-RP1. CD275 is a ligand for ICOS (CD278) that is expressed on activated T cells. CD275 is expressed by antigen presenting cells including B cells, macrophages, and dendritic cells. The interaction of CD275 with its receptor ICOS plays an important role in T cell costimulation.

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	Assay dependent
Inhibition assay	Assay dependent
Neutralizing /Inhibition	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	10mM NaH ₂ PO ₄ , 150mM NaCl
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Purification	Purified by affinity chromatography From tissue culture supernatant
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).