

Integrin alpha 4 + beta 7 antibody [DATK32] (PE)

Cat. No. GTX01493-08

Host	Rat
Clonality	Monoclonal
Isotype	IgG2a
Applications	FCM
Reactivity	Mouse

References (3)

Package

100 µg

PRODUCT

Summary

The DATK32 antibody is specific for a combinatorial epitope of the mouse integrin α 4/ β 7 (LPAM-1) heterodimer. LPAM-1 is expressed on most mature lymphocytes and sub-populations of thymocytes and bone marrow cells. It is involved in lymphocyte adhesion and migration and is a receptor for several extracellular matrix proteins including VCAM-1, fibronectin, and MAdCAM-1. DATK32 is reported to induce LPAM-1 dependent lymphocyte aggregation and block adhesion to VCAM-1, fibronectin, and MAdCAM-1.

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution

Recommended dilution

FCM

Assay dependent

Not tested in other applications.

Properties

Form Liquid

Buffer 10mM NaH₂PO₄, 150mM NaCl, 0.1% Gelatin

Preservative 0.09% Sodium azide

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.

Concentration 0.2 mg/ml (Please refer to the vial label for the specific concentration.)

Purification Purified by affinity chromatography
From tissue culture supernatant

Conjugation Phycoerythrin (PE) [Wavelength](#)

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

**GTX01493-08 FCM Image**

FACS analysis of mouse C57Bl/6 bone marrow cells using GTX01493-08 Integrin alpha 4 + beta 7 antibody [DATK32] (PE).

Right panel : co-stained with Integrin alpha 4 + beta 7 antibody [DATK32] (PE) and Mouse CD45R antibody (APC)

Left panel : co-stained with isotype control and Mouse CD45R antibody (APC)
antibody amount : 0.5 μ g (5 μ l)



For full product information, images and publications, please visit our [website](#).