

Armenian Hamster IgG Isotype Control [PIP] (APC)

Cat. No. GTX01517-07

Host Armenian Hamster**Clonality** Monoclonal**Isotype** IgG**Applications** FCM

Package

100 µg

PRODUCT

Summary

The PIP immunoglobulin is useful as an isotype-matched control. The PIP immunoglobulin is specific for Glutathione S-transferase and can be used as an isotype control for Armenian Hamster IgG antibodies.

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution

FCM

Recommended dilution

Assay dependent

Not tested in other applications.

Properties

Form Liquid**Buffer** 10mM NaH₂PO₄, 150mM NaCl, 0.1% Gelatin**Preservative** 0.09% Sodium azide**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.**Concentration** 0.2 mg/ml (Please refer to the vial label for the specific concentration.)**Purification** Purified by affinity chromatography
From tissue culture supernatant**Conjugation** Allophycocyanin (APC) [Wavelength](#)

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).