

## Armenian Hamster IgG Isotype Control [PIP] (PE)

Cat. No. GTX01517-08

**Host** Armenian Hamster**Clonality** Monoclonal**Isotype** IgG**Applications** FCM

Package

100 µg

## PRODUCT

## Summary

The PIP immunoglobulin is useful as an isotype-matched control. The PIP immunoglobulin is specific for Glutathione S-transferase and can be used as an isotype control for Armenian Hamster IgG antibodies.

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

## Suggested dilution

## Recommended dilution

FCM

Assay dependent

Not tested in other applications.

## Properties

**Form** Liquid**Buffer** 10mM NaH<sub>2</sub>PO<sub>4</sub>, 150mM NaCl, 0.1% Gelatin**Preservative** 0.09% Sodium azide**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.**Concentration** 0.2 mg/ml (Please refer to the vial label for the specific concentration.)**Purification** Purified by affinity chromatography  
From tissue culture supernatant**Conjugation** Phycoerythrin (PE) [Wavelength](#)

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).