

## SARS-CoV-2 (COVID-19) 3CLpro (nsp5) protein, His and Avi tag

**Cat. No. GTX01557-pro****Species** SARS Coronavirus 2**Package**  
100 µg


## Properties

**Form** Liquid**Buffer** Filter-sterilized 20mM Tris, 250mM NaCl, 0.5% Triton X-100**Preservative** No preservatives**Storage** Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.**Region/Sequence** N-terminal 6×His and Avi tagged; Ser1-Gln306 of 2019-ncov 3C-like Proteinase (Accession #YP\_009725301.1)**Expression System** E. coli**Purity** > 95% by SDS-PAGE**Endotoxin** < 1.0 EU per µg (determined by the LAL method)**Conjugation** Unconjugated**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES

**GTX01557-pro Image**

SDS-PAGE of GTX01557-pro SARS-CoV-2 (COVID-19) 3C-like Proteinase protein, His and Avi tag.



For full product information, images and publications, please visit our [website](https://www.genetex.com).