

# SARS-CoV-2 (COVID-19) Envelope Protein, GST and His Tag

**Cat. No. GTX01565-pro**
**Species** SARS Coronavirus 2

**Package**  
100 µg

## PROPERTIES

**Form** Lyophilized powder

**Buffer** Batch dependent (please contact us for details)

**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, we recommend adding a carrier protein (0.1% HSA or BSA), then aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

**Region/Sequence** Coronavirus 2019 CoV-2 Envelope, Wuhan-Hu-1 strain, amino acids 1-75 fused to GST-tag at N-terminal & His-tag at C-terminal.

**Expression System** HEK293 cells

**Purification** Purified by Metal-Afinity chromatographic technique

**Purity** > 90% by SDS-PAGE

**Conjugation** Unconjugated

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


For full product information, images and publications, please visit our [website](#).