

## beta Catenin antibody [17C2]

Cat. No. GTX01865

|              |            |
|--------------|------------|
| Host         | Mouse      |
| Clonality    | Monoclonal |
| Isotype      | IgG2a      |
| Applications | IHC-P      |
| Reactivity   | Human      |

**Package**  
500 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| IHC-P              | 1:150                |

**Note : The recommended fixative is 10% neutral-buffered formalin for paraffin-embedded tissue sections. Heat Induced Epitope Retrieval (HIER) was performed in Epitope Retrieval Solution pH 9.**

Not tested in other applications.

## Properties

|              |                                                                                                                          |
|--------------|--------------------------------------------------------------------------------------------------------------------------|
| Form         | Liquid                                                                                                                   |
| Buffer       | Tissue culture supernatant                                                                                               |
| Preservative | 0.09% Sodium azide                                                                                                       |
| Storage      | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.                   |
| Immunogen    | Prokaryotic recombinant protein corresponding to a 160 amino acid region of the C-terminus of the beta-catenin molecule. |
| Purification | Unpurified                                                                                                               |
| Conjugation  | Unconjugated                                                                                                             |


## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

**DATA IMAGES**

**GTX01865 IHC-P Image**

IHC-P analysis of human colon polyp using GTX01865 beta Catenin antibody [17C2]. Note the abnormal translocation of the protein to the nucleus.



For full product information, images and publications, please visit our [website](#).