

CD79b antibody [JS01]

Cat. No. GTX01875

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	IHC-P
Reactivity	Human

Package
500 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:40

Note : High temperature antigen unmasking technique. Standard ABC technique.

Not tested in other applications.

Properties

Form	Liquid
Buffer	Tissue culture supernatant
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Immunogen	A recombinant protein containing 48 amino acids of the cytoplasmic domain of the CD79b molecule.
Purification	Unpurified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX01875 IHC-P Image

IHC-P analysis of human tonsil tissue using GTX01875 CD79b antibody [JS01]. Note intense membrane staining of B cells.



For full product information, images and publications, please visit our [website](#).