

## Calcitonin antibody [CL1948]

Cat. No. GTX01932

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	IHC-P
Reactivity	Human

Package  
500 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:200

**Note : Enzyme Induced Epitope Retrieval (EIER) was performed in Enzyme Proteinase K.**

Not tested in other applications.

## Properties

Form	Liquid
Buffer	Tissue culture supernatant
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Immunogen	Prokaryotic recombinant fusion protein containing 32 amino acids corresponding to the mature human calcitonin molecule.
Purification	Unpurified
Conjugation	Unconjugated

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTx01932 IHC-P Image

IHC-P analysis of human medullary thyroid carcinoma tissue using GTx01932 Calcitonin antibody [CL1948].



For full product information, images and publications, please visit our [website](#).