

beta Sarcoglycan antibody [β Sarc1/5B1]**Cat. No. GTX01942**

| | |
|---------------------|------------|
| Host | Mouse |
| Clonality | Monoclonal |
| Isotype | IgG1 |
| Applications | IHC-Fr |
| Reactivity | Human |

Package
500 μ l

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| IHC-Fr | 1:50–1:200 |

Note : Freeze specimen tissue blocks in isopentane chilled in liquid nitrogen.

Not tested in other applications.

Properties

| | |
|---------------------|--|
| Form | Liquid |
| Buffer | Tissue culture supernatant |
| Preservative | 0.09% Sodium azide |
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. |
| Immunogen | Fusion protein RBSG-NT of the human beta-sarcoglycan sequence. |
| Purification | Unpurified |
| Conjugation | Unconjugated |

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).