

delta Sarcoglycan antibody [δ Sarc3/12C1]**Cat. No. GTX01943**

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	IHC-Fr
Reactivity	Human

Package
500 μ l

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-Fr	1:25–1:50

Note : Freeze specimen tissue blocks in isopentane chilled in liquid nitrogen.

Not tested in other applications.

Product Note We do not recommend use of this product for Mouse,Rat,Rabbit,Dog,Hamster,Chicken,Pig samples.

Properties

Form	Liquid
Buffer	Tissue culture supernatant
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Synthetic peptide containing amino acids 1-19 at the N-terminus of the human delta-sarcoglycan sequence (Jung D. et al.,1966).
Purification	Unpurified
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).