

# delta Sarcoglycan antibody [ $\delta$ Sarc3/12C1]

**Cat. No. GTX01943**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2a
<b>Applications</b>	IHC-Fr
<b>Reactivity</b>	Human

**Package**  
500  $\mu$ l

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-Fr	1:25–1:50

**Note : Freeze specimen tissue blocks in isopentane chilled in liquid nitrogen.**

Not tested in other applications.

### Product Note

We do not recommend use of this product for Mouse,Rat,Rabbit,Dog,Hamster,Chicken,Pig samples.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	Tissue culture supernatant
<b>Preservative</b>	0.09% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Immunogen</b>	Synthetic peptide containing amino acids 1-19 at the N-terminus of the human delta-sarcoglycan sequence (Jung D. et al.,1966).
<b>Purification</b>	Unpurified
<b>Conjugation</b>	Unconjugated

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).