

gamma Sarcoglycan antibody [35DAG/21B5]

Cat. No. GTX01944

| | |
|--------------------|---------------|
| Host | Mouse |
| Clonality | Monoclonal |
| Isotype | IgG2b |
| Application | IHC-Fr |
| Reactivity | Human, Rabbit |

Package
500 µl

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| IHC-Fr | 1:100 |

Note : Freeze specimen tissue blocks in isopentane chilled in liquid nitrogen.

Not tested in other applications.

Calculated MW 32 kDa. ([Note](#))

PROPERTIES

| | |
|---------------------|--|
| Form | Liquid |
| Buffer | Tissue culture supernatant |
| Preservative | 0.09% Sodium azide |
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Immunogen | Synthetic peptide containing amino acids 167–178 of the rabbit gamma-sarcoglycan sequence (Noguchi et al., 1995). |
| Purification | Unpurified |
| Conjugation | Unconjugated |

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).