

Human MMP7 protein (Active)

Cat. No. GTX02553-pro

Applications	Functional Assay
Species	Human

Package
20 µg

Applications

Application Note

The specific activity was found to be 6,200 U/mg. One unit is defined as the digestion of 1 µg Azocoll/min at 37°C.

Properties

Form	Liquid
Buffer	10mM HEPES, 5mM CaCl ₂ , 150mM NaCl
Preservative	No preservatives
Storage	Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.
Expression System	E. coli
Purification	Purified by proprietary chromatographic techniques
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).