

CDC4 antibody

Cat. No. GTX02557

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

References (1)
Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	2.5 µg/ml

Not tested in other applications.

Calculated MW 80 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 2% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5-1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of Human FBXW7.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX02557 IHC-P Image

IHC-P analysis of human pineal tissue using GTX02557 CDC4 antibody. Cytoplasmic in cell bodies and processes of pinealocytes.

Dilution : 1:100



GTX02557 IHC-P Image

IHC-P analysis of human kidney tissue using GTX02557 CDC4 antibody.


GTX02557 WB Image

WB analysis of human fetal heart tissue lysate using GTX02557 CDC4 antibody.

Dilution : 1 µg/ml


GTX02557 WB Image

WB analysis of Jurkat cell lysate using GTX02557 CDC4 antibody.

Dilution : 0.2-1 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).