

ER81 antibody

Cat. No. GTX02562

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent

Not tested in other applications.

Calculated MW 55 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 2% Sucrose
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 - 1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human ETV1.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES



GTX02562 IHC-P Image

IHC-P analysis of human cortex tissue using GTX02562 ER81 antibody.


Dilution : 5 µg/ml



GTX02562 WB Image

WB analysis of MDA-MB-435s cell lysate using GTX02562 ER81 antibody.

Dilution : 1 µg/ml



GTX02562 WB Image

WB analysis of ER81-transfected HEK293 cell lysate using GTX02562 ER81 antibody.

Lane 1 : 100 µg untransfected HEK293 lysate

Lane 2 : 33 µg ER81 transfected HEK293 lysate

Dilution : 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).