

Mouse Anti-Armenian Hamster IgG1 antibody [SB139a] (HRP)

Cat. No. GTX02569

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, FCM, ELISA
Reactivity	Armenian Hamster, Golden Syrian Hamster

Package
500 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
FCM	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Product Note Armenian/Syrian hamster IgG1

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Armenian Hamster IgG
Purification	Purified IgG1
Conjugation	Horseradish peroxidase(HRP)


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX02569 ELISA Image

ELISA analysis of purified Armenian hamster IgG1, IgG2, and IgG3 using serially diluted GTX02569 Mouse Anti-Armenian Hamster IgG1 antibody [SB139a] (HRP).



For full product information, images and publications, please visit our [website](https://www.genetex.com).