

Borrelia burgdorferi VlsE1 protein, His tag

Cat. No. GTX02573-pro

Species	Borrelia burgdorferi
----------------	----------------------

Package	20 µg
----------------	-------

Properties

Form	Liquid
-------------	--------

Buffer	20mM HEPES, 200mM NaCl, 20% Glycerol
---------------	--------------------------------------

Preservative	No preservative
---------------------	-----------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
----------------	--

Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
----------------------	--

Region/Sequence	N-terminal 10x His tagged Borrelia VlsE1 protein
------------------------	--

Expression System	E. coli
--------------------------	---------

Purification	Purified by proprietary chromatographic techniques
---------------------	--

Purity	> 95% by SDS-PAGE
---------------	-------------------

Conjugation	Unconjugated
--------------------	--------------

Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>
-------------	---



For full product information, images and publications, please visit our [website](#).