

CD74 antibody [CLIP/3127R]

Cat. No. GTX02612

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IHC-P, FCM
Reactivity	Human, Mouse, Baboon

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1-2 µg/ml
IHC-P	1-2 µg/ml
FCM	1-2 µg/million cells

Note : Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris buffer with 1mM EDTA (pH 9.0) for 45 min at 95°C followed by cooling at RT for 20 minutes.

Not tested in other applications.

Calculated MW	34 kDa. (Note)
Observed MW (kDa)	~40 kDa. The observed M.W. is based on the following publication. PMID: 29490700
Product Note	We do not recommend use of this product for Rat samples.

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	200 µg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	SU-DHL-4 lymphoma cells
Purification	Protein A/G purified
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX02612 IHC-P Image

IHC-P analysis of human tonsil tissue section using GTX02612 CD74 antibody [CLIP/3127R].




GTX02612 FCM Image

FACS analysis of Raji cells using GTX02612 CD74 antibody [CLIP/3127R].

Blue : Primary antibody

Red : Isotype control



GTX02612 WB Image

WB analysis of 293T whole cell lysate using GTX02612 CD74 antibody [CLIP/3127R].



For full product information, images and publications, please visit our [website](#).