

# IL10 antibody [IL10/2651R]

**Cat. No. GTX02650**

|                     |            |
|---------------------|------------|
| <b>Host</b>         | Rabbit     |
| <b>Clonality</b>    | Monoclonal |
| <b>Isotype</b>      | IgG        |
| <b>Applications</b> | IHC-P      |
| <b>Reactivity</b>   | Human      |

**Package**  
100 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution                                                                                                                                                                        | Recommended dilution |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------|
| IHC-P                                                                                                                                                                                     | 1-2 µg/ml            |
| <b>Note : Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris buffer with 1mM EDTA (pH 9.0) for 45 min at 95°C followed by cooling at RT for 20 minutes.</b> |                      |

Not tested in other applications.

## Properties

|                      |                                                                                         |
|----------------------|-----------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                  |
| <b>Buffer</b>        | PBS, 0.05% BSA                                                                          |
| <b>Preservative</b>  | 0.05% Sodium azide                                                                      |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. |
| <b>Concentration</b> | 200 µg/ml (Please refer to the vial label for the specific concentration.)              |
| <b>Immunogen</b>     | Recombinant full-length human Interleukin 10 (IL10) protein                             |
| <b>Purification</b>  | Protein A purified                                                                      |
| <b>Conjugation</b>   | Unconjugated                                                                            |

### Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTx02650 IHC-P Image

IHC-P analysis of human endocervical carcinoma section using GTx02650 IL10 antibody [IL10/2651R].



### GTx02650 IHC-P Image

IHC-P analysis of human prostate carcinoma section using GTx02650 IL10 antibody [IL10/2651R].



For full product information, images and publications, please visit our [website](#).