

## Human FGFR2 alpha IIIb (extracellular region) protein, His tag

**Cat. No. GTX02783-pro****Species** Human**Package**

50 µg

## Applications

**Observed MW (kDa)** ~80 kDa.


## Properties

**Form** Liquid**Buffer** PBS, 20% Glycerol**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (at least 1 month), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** batch dependent (Please refer to the vial label for the specific concentration.)**Region/Sequence** C-terminal 10 x His tagged ; 22-378 aa (Uniprot #P21802-3)**Expression System** HEK293 cells**Purity** > 95% by SDS-PAGE**Endotoxin** < 1 EU/µg (Determined by LAL assay)**Conjugation** Unconjugated**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX02783-pro Image

SDS-PAGE of GTX02783-pro Human FGFR2 alpha IIIb (extracellular region) protein.



For full product information, images and publications, please visit our [website](https://www.genetex.com).

Date 2026 / 01 / 31 Page 2 of 2