

## Human PLA2R1 (extracellular region) protein, His tag

Cat. No. GTX02788-pro

Species	Human
---------	-------

Package	50 µg
---------	-------

## Applications

Observed MW (kDa)	~200 kDa.
-------------------	-----------

## Properties

Form	Liquid
------	--------

Buffer	PBS, 20% Glycerol
--------	-------------------

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (at least 1 month), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Concentration	batch dependent (Please refer to the vial label for the specific concentration.)
---------------	----------------------------------------------------------------------------------

Region/Sequence	C-terminal 10 x His tagged ; 21-1397 aa (Uniprot #Q13018)
-----------------	-----------------------------------------------------------

Expression System	HEK293 cells
-------------------	--------------

Purity	> 95% by SDS-PAGE
--------	-------------------

Endotoxin	< 1 EU/µg (Determined by LAL assay)
-----------	-------------------------------------

Conjugation	Unconjugated
-------------	--------------


Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
------	---------------------------------------------------------------------------------------------------------------------------------------------------

Note	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTx02788-pro Image

SDS-PAGE of GTx02788-pro Human PLA2R1 (extracellular region) protein.



For full product information, images and publications, please visit our [website](https://www.genetex.com).