

Human PLA2R1 (extracellular region) protein, His tag

Cat. No. GTX02791-pro

Species Human

Package
50 µg

Applications

Observed MW (kDa) ~24 kDa.

Properties

Form Liquid

Buffer PBS, 20% Glycerol

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (at least 1 month), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration batch dependent (Please refer to the vial label for the specific concentration.)

Region/Sequence C-terminal 10 x His tagged ; 21-173 aa (Uniprot #Q13018)

Expression System HEK293 cells

Purity > 95% by SDS-PAGE

Endotoxin < 1 EU/µg (Determined by LAL assay)

Conjugation Unconjugated


Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx02791-pro Image

SDS-PAGE of GTx02791-pro Human PLA2R1 (extracellular region) protein.



For full product information, images and publications, please visit our [website](#).