

Human VAP1 (extracellular region) protein, His tag

Cat. No. GTX02794-pro

Species	Human
---------	-------

Package	50 µg
---------	-------

Applications

Observed MW (kDa)	~100 kDa.
-------------------	-----------

Properties

Form	Liquid
------	--------

Buffer	PBS, 20% Glycerol
--------	-------------------

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (at least 1 month), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	---

Concentration	batch dependent (Please refer to the vial label for the specific concentration.)
---------------	--

Region/Sequence	N-terminal 10 x His tagged ; 29-763 aa (Uniprot #Q16853)
-----------------	--

Expression System	HEK293 cells
-------------------	--------------

Purity	> 95% by SDS-PAGE
--------	-------------------

Endotoxin	< 1 EU/µg (Determined by LAL assay)
-----------	-------------------------------------

Conjugation	Unconjugated
-------------	--------------

Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
------	---

Note	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
------	---



For full product information, images and publications, please visit our [website](#).