

HAI-2 antibody [XY9]

Cat. No. GTX02801

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Human

References (1)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1-5 µg/ml
ICC/IF	Assay dependent

Not tested in other applications.

Calculated MW 28 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant human protein expressed in mammalian cells
Purification	Purified by protein A
Purity	>90% determined by SDS-PAGE
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX02801 WB Image

Untreated (–) and treated (+) MCF-7 whole cell extracts (30 µg) were separated by 12% SDS-PAGE, and the membrane was blotted with HAI-2 antibody [XY9] (GTX02801) diluted at 1:1000. The HRP-conjugated anti-mouse IgG antibody (GTX213111-01) was used to detect the primary antibody.

Double bands were observed at 25~45 kDa.

Reference: PMID: 25786220



GTX02801 WB Image

Unboiled and boiled various whole cell extracts (30 µg) were separated by 12% SDS-PAGE, and the membrane was blotted with HAI-2 antibody [XY9] (GTX02801) diluted at 1:1000. The HRP-conjugated anti-mouse IgG antibody (GTX213111-01) was used to detect the primary antibody.

Double bands were observed at 25~45 kDa.

Reference: PMID: 25786220



For full product information, images and publications, please visit our [website](https://www.genetex.com).