

DDX17 antibody [GT1254]

Cat. No. GTX02851

Host	Rabbit
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
IHC-P	1:50 - 1:200

Not tested in other applications.

Calculated MW	80 kDa. (Note)
Observed MW (kDa)	72 82 kDa. The observed M.W. is based on the following publication. PMID: 31653828

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, we recommend adding 50% glycerol, then aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.54 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthesized peptide derived from human DDX17
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx02851 WB Image


Various whole cell extracts (30 µg) were separated by 7.5% SDS-PAGE, and the membrane was blotted with DDX17 antibody [GT1254] (GTx02851) diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (GTx213110-01) was used to detect the primary antibody.
The observed M.W. is based on the publication: PMID: 31653828


GTx02851 IHC-P Image

IHC-P analysis of mouse spinal cord section using GTx02851 DDX17 antibody [GT1254].
Dilution : 1:100


GTx02851 IHC-P Image

IHC-P analysis of rat testis tissue section using GTx02851 DDX17 antibody [GT1254].
Dilution : 1:100


GTx02851 WB Image

WB analysis of various samples using GTx02851 DDX17 antibody [GT1254].
Dilution : 1:1000
Loading : 25µg



For full product information, images and publications, please visit our [website](https://www.genetex.com).