

Human LBP protein, His tag

Cat. No. GTX02867-pro


Species	Human
---------	-------

Package	50 µg
---------	-------

Properties

Form	Liquid
Buffer	PBS, 10% Glycerol
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1 week), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.25 mg/ml (Please refer to the vial label for the specific concentration.)
Region/Sequence	His-Tagged; 26-481aa (P18428, NP_004130.2), signal peptide removed
Expression System	HEK293 cells
Purity	> 95% by SDS-PAGE
Endotoxin	< 1.0 EU/µg (determined by LAL method)
Conjugation	Unconjugated
Note	For <i>In vitro</i> laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

DATA IMAGES



GTX02867-pro Image

SDS-PAGE analysis of 3µg GTX02867-pro Human LBP protein, His tag under reducing condition.

For full product information, images and publications, please visit our [website](#).