

Apolipoprotein E antibody [E276]

Cat. No. GTX02904

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IP, ELISA, Sandwich ELISA
Reactivity	Human

Package
250 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IP	Assay dependent
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture: GTX02904 / GTX02906, Detection: GTX02905-02.

Not tested in other applications.

Calculated MW 36 kDa. ([Note](#))**Product Note** This antibody is able to detect native and recombinant human apolipoprotein E (isoforms apoE2, apoE3 and apoE4).

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified human apolipoprotein E
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated


For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx02904 ELISA Image

Sandwich ELISA analysis of human apoE protein using GTx02904 Apolipoprotein E antibody [E276] as coating antibody and GTx02905-02 Apolipoprotein E antibody [E887] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).