

# Apolipoprotein E antibody [mE29] (Biotin)

**Cat. No. GTX02908-02**

<b>Host</b>	Rat
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG2a
<b>Applications</b>	WB, ELISA, Sandwich ELISA
<b>Reactivity</b>	Mouse

**Package**  
250 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

**Note : Capture: GTX02907, Detection: GTX02908-02.**

Not tested in other applications.

**Calculated MW** 36 kDa. ( [Note](#) )

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.02% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Recombinant mouse apoE
<b>Purification</b>	Protein G purified From tissue culture supernatant
<b>Conjugation</b>	Biotin Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

## Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES



### GTx02908-02 ELISA Image

Sandwich ELISA analysis of mouse apoE protein using GTx02907 Apolipoprotein E antibody [mE1] as coating antibody and GTx02908-02 Apolipoprotein E antibody [mE29] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).