

## Interferon gamma antibody [MT17.1]

Cat. No. GTX02935

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Application	ELISA, Sandwich ELISA
Reactivity	Sheep, Bovine

Package  
250 µg

## APPLICATION

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

**Note : Capture: GTX02935, Detection: GTX02936-02.**

Not tested in other applications.

Calculated MW	19 kDa. ( <a href="#">Note</a> )
Product Note	This antibody is able to detect native and recombinant bovine IFN-γ.

## PROPERTIES

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Native and recombinant bovine IFN-γ
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

Date 2024 / 05 / 17 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

**Note**

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

**DATA IMAGES**

**GTx02935 ELISA Image**

Sandwich ELISA analysis of bovine IFN-gamma protein using GTx02935 Interferon gamma antibody [MT17.1] as coating antibody and GTx02936-02 Interferon gamma antibody [MT307] (Biotin) as detecting antibody.

Substrate : pNPP



For full product information, images and publications, please visit our [website](https://www.genetex.com).