

Interferon gamma antibody [MT327]

Cat. No. GTX02938

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ELISA, Sandwich ELISA
Reactivity	Rabbit

Package
250 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture: GTX02938, Detection: GTX02937-02.

Not tested in other applications.

Product Note This antibody is able to detect native and recombinant rabbit IFN- γ .

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant rabbit IFN- γ
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx02938 ELISA Image

Sandwich ELISA analysis of rabbit IFN-gamma protein using GTx02938 Interferon gamma antibody [MT327] as coating antibody and GTx02937-02 Interferon gamma antibody [MT318] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).